Citieffe srl a socio unico
Via Armaroli, 21
40012 Calderara di Reno (Bologna) - Italy
Tel. +39 051 721850 - Fax +39 051 721870
info@citieffe.com - www.citieffe.com

ITALIANO

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ENGLISH

INSTRUCTIONS FOR USE OF X-FRAME EXTERNAL FIXATORS

Please, read before use.

This package insert does not include all of the information necessary for selection and use of a device. Please see full labelling for all necessary information.

Federal law (U.S.A.) restricts this device to be sale by or on the order of a licensed physician.

1. Description and intended use

X-Frame external fixators are single-use medical devices to be used as an aid in the treatment of bone fracture repair and healing. The range includes medical devices from Class I, sterile

compliant with EU 2017/405.

The system is composed of clamps, rods and accessories designed to build an external fixator construct. For further information regarding models, dimensions and sizes see the Citieffe catalogues.

2. Materials

The materials used for devices manufacturing are:

- aluminium alloy EN AW-7075 (ISO 573-3 standard);
- stainless steels: AISI316L (ASTM F-138, ISO 5832-1 standards), AISI303 (EN 10088-3 standard);
- Ti6Al4V El titanium alloy (ASTM F-136, ISO 5832-3 standards);
- PEEK CA-30 radiolucent technopolymer;
- carbon-fiber and epoxy-resin composite;
- silicone.

3. Sterilization

The devices are provided sterile and have been exposed to validated ethylene oxide (EO) sterilization process.

Sterile devices shall **not** be sterilized or resterilized by the user before use.

The manufacturer and the distributor accept no liability for sterilization of the devices performed by the user.

4. General safety and performance notes

- For the operative technique see the surgical technique manual of each product.
- The devices are designed according to state-of-the-art knowledge of implant technology and the applicable rules of engineering skills. However, the safety and functionality of devices can only be guaranteed if surgeons observe and comply with the notes in this instruction for use and the surgical technique.
- The sole purpose of devices is to support the healing process. They are no replacement for intact tissue and bone material. Due to the anatomy of human bones, implants are subject to limits with regard to their size and thickness.
- The tolerance, the resistance and the duration of the system depend on the material, on the correct use and assembling of the components and on the use of the dedicated instrument set.

5. Storage, handling and treatment of devices

- Always check the expiration date (shelf life) and the indicator of sterility on the sterile packaging before opening. If the expiration date has past or if the indicator of sterility indicates non-sterile, the devices should be considered non-sterile and unusable.
- The sterile packaging should be checked for defects in the sterile barrier before opening. If any defects are noticed, the devices should be considered non-sterile and unusable.
- Open the sterile packaging device only during surgery, in the sterile field and immediately before use. Follow the hospital protocol in order to maintain the sterility of the device during handling. Remove all protective packaging before implantation.
- Visually inspect the devices before use. Any devices with damage, scratches, flaws, corrosion, discoloration, residue or debris should be discarded.
- Use extreme care in the handling of devices. Cutting, bending or scratching the surface of metal components can significantly reduce the strength and fatigue resistance of the device.
- Do not machine or modify devices in any way, unless required by the surgical technique and allowed by the manufacturer. In case of doubt, request a written recommendation from the manufacturer.
- Store products packed in their original packaging.
- Users should always wear appropriate personal protective equipment when handling the devices.
- Storage areas for devices should be clean from dust, away from areas of humidity, away from sunlight and at room temperature. This recommendation is equally valid for the transport and packaging of the devices.

6. Intended users

- The devices shall only be handled and used by qualified and trained healthcare practitioner (e.g. surgeons who are experienced in orthopedic surgery, operating room staff, and individuals involved in preparation of the device).
- Surgery shall only be performed by operating surgeons who are experienced in orthopedic surgery and who have acquired the appropriate qualifications, are aware of general risks of orthopedic surgery and are familiar with the product-specific surgical procedures.
- The surgeon's education, training and professional judgment must be relied upon to choose the most appropriate implant type and size and treatment option. No mandatory training on the product is required before using X-Frame external fixators.
- Citieffe cannot be held liable for complications arising from incorrect indication, surgical techniques or asepis; this is the sole responsibility of the surgeon.

7. Target group of patients

The intended target group for the treatment of temporary fracture fixation with external fixators consists of all patients who have a clinical condition that can be attributed to one of the indications for use of the X-Frame external fixators.

From preclinical evaluation, state of the art, and post market surveillance data, there are no known age-related limitations of use, so special populations such as paediatric (0-17) and elderly (>65) are also included in the target patient groups.

However, it is not recommended to use the device on patients who have one or more of the contraindications listed in these instructions for use. Moreover, further precautions, especially in vulnerable populations (e.g. pregnant and breast-feeding populations) or suspected contraindicated population, should be taken as indicated in section "Suggestions" 11, keeping in mind also the risks outlined in section 12.

8. Indications

X-Frame external fixators are indicated for the treatment of open or closed fractures in long bones, pelvic fractures, infected and aseptic non-unions, arthrodesis.

9. Contraindications

- In progress infection or infection risk due to present or previous pathologies in the area to be treated.
- Conditions which tend to pre-empt the patient's ability or willingness to restrict activities or follow direction during treatment.
- Bone stock compromised by disease, infection or prior implantation that cannot provide adequate support and/or fixation of the devices.
- Any mental or neuromuscular disorder which could create an unacceptable risk of fixation failure or complications in postoperative care.
- Any drug addiction and/or tendency to abuse drugs, medicines or alcohol.
- Other medical or surgical conditions which would preclude the potential benefit of surgery.

10. Warnings

- Proper understanding of the devices and surgical technique are essential. Physicians are strongly encouraged to obtain instruction from experienced clinicians and to observe surgical applications of the apparatus prior to initial use.
- Repeat use of a single-use device is strictly forbidden. Each device used once shall be disposed of according to the hospital procedures. This is the same even where it appears to be intact; the device may have small faults or internal stresses that if the item is used may lead to fatigue failure.
- It is strictly prohibited to carry out any modification whatsoever on devices. If this recommendation is not followed, Citieffe disclaims any liability for any subsequent consequences.
- It is strictly forbidden to recondition or resterilize sterile devices taken out of the original packaging. Any devices that have been taken out of its sterile package but have not been used or have become soiled should be handled according to hospital protocol.
- Do not use damaged, scratched, improperly treated devices or devices handled by unauthorized persons.
- Do not use any kind of disinfectants and powder chemicals on the surface of the fixator bodies, components and accessories.

humans, implants are still subjected to limitations in cuanto a tamaño y espesor.

Mixing of devices from different suppliers is not recommended for reasons of metallurgy, mechanics and design unless otherwise specified (see the surgical technique). Citieffe disclaims any liability in case of implants from different sources being mixed.

11. Suggestions

- A preliminary review of the special surgical instruments is recommended.
- Use only original Citieffe devices and accessories.
- It's extremely important to select the most suitable type of fixator and bone screws to be used depending on bone.
- The components of the fixator may be of radiolucent technopolymer and carbon fiber and epoxy-resin composite should not be sterilized in a steam autoclave.
- The medical staff should inform the patient of the appropriate pin site care required, and the right adjustments to be periodically done on the device and its accessories.
- At the end of the treatment, once the devices have been removed, they should be disposed of as medical special waste.

12. Residual risks, adverse events and possible undesirable side effects

- Device bending or breakage due to an excessive intra or post-operative stress.
- Fixator clamps loosening.
- Mobilization and loss of bone fixation.
- Limb length inequality caused by devices poor positioning or choice.

There are risks risks associated with anesthesia and surgery and residual risks related to the failure of the device or complications. The latter risks may be caused by improper use of the device with respect to this manual and the intended surgical technique, use in the presence of contraindications, incorrect device selection and placement, reuse of a disposable device. The surgeon and the patient must be aware of these risks and possible need of an early revision surgery.

Any incidents associated with the use of the product, including those other than those mentioned above, shall be reported without undue delay to the manufacturer and the national authority.

Manufacturer's contact for reporting:

- support@citieffe.com

13. Patient information

- The patient must be informed of the residual risks, adverse events and possible undesirable side effects associated with the use of the medical device, as well as the implications for carrying capacity, mobility and general living conditions.
- The patient must be informed of the possibility of a refracture occurring even when the device is in place or immediately after its removal; the possibility is less later, when the space left free by the implant is filled.
- The patient must be made aware of and properly instructed by the medical staff regarding medication in the areas close to the bone screw insertion point, postoperative activities and the need to attend postoperative check-ups.

14. MRI Safety Information

Non-clinical testing has demonstrated that X-Frame external fixator is MR Conditional. It can be scanned safely under the following conditions:

- Static magnetic field: 1.5T and 3T;
 - Maximum spatial field gradient: 1100 Gauss/cm (11 T/m) or less;
 - Maximum value of the specific absorption rate (SAR) reported by the MR system as a whole-body average: 2W/Kg;
 - X-Frame external fixator must be outside the MR scanner bore.
- Under the scan conditions defined above the X-Frame external fixators do not produce a temperature rise at 1.5T and 3T after 15 minutes of continuous scanning.

Non-clinical testing has demonstrated that X-Frame external fixators scanned with static magnetic field of 1.5T and 3T showed:

- No displacement force and torque effects;
 - Image artifacts up to 70 mm from the screws using Spin Echo sequence and up to 39 mm from the screws using Gradient Echo sequence.
- For further information on tested components and clinical configurations, please refer to the surgical technique of the medical device.

For detailed information regarding the safety and compatibility of the devices in the MR (Magnetic Resonance) environment, the Physician should contact Citieffe.

15. Traceability

To ensure the traceability of the devices, the user should attach the labels provided in the packaging to the patient's clinical records.

16. Symbols

Symbols in the label are harmonized and in accordance with EN ISO 15223-1.

Non-harmonized symbols:

R_x Only "Caution: Federal law (USA) restricts this device to sale by or on the order of a physician"

The symbols glossary is available at this address <https://www.citieffe.com/symbolsglossary>

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